PERFORMANCE EVALUATION: OXYGEN CONCENTRATORS

Date: Lab	Clinical		_ Agency		
Lab: Pass F		Clinical: Pass			
Student name		Instructor name .			
No. of times observ	ed in clinical				
No. of times practic	ed in clinical				
PASSING CRITERIA: Obtain 90% or better on the pro- the evaluation is terminated. Pro- performance receives a failing g		ocedure must be per	ed by * must re formed within	ceive at lea designated	st 1 point, or time, or the
SCORING:	2 points — Task performed satist 1 point — Task performed satist 0 points — Task performed inco NA — Task not applicable to the	actorily with self-initi rrectly or with promp	iated correction oting required.	1.	
TASKS:			PEER	LAB	CLINICAL
1. Places the	e concentrator correctly at the poi	nt of use		_	
* a. Away from walls or draperies					
* b. Away from heaters, heat registers, or radiators		adiators			
 Connects the power cord, ensuring that the minimal electrical load 		e circuit has			
* 3. Checks th	ne gross particle filter and cleans a	as needed			
* 4. Connects a humidifier as required					Ц
* 5. Connects the oxygen delivery device					Ц
* 6. Adds up to 50 feet of extension tubing as needed		needed			
* 7. Turns on the power switch					
* 8. Sets the liter flow rate				Ц	
* 9. Documents procedure in the patient record					
SCORE: Peer	points of possible 2	0;%			
Lab	points of possible 2	0;%			
Clinica	points of possible 2	0;%			
TIME: ou	t of possible 30 minutes				
STUDENT SIGNATU	JRES	INSTRUCTOR SI	GNATURES		
PEER:		LAB:			
STUDENT:		CLINICAL:			